
 Keim GmbH Rheinstraße 1-5, 63225 Langen Standort: Otto-Hahn-Straße 41, 63456 Hanau, Tel. +49 6181-9886-10	Job No.: me511042	Created at: 15.10.2019	Operator: bka
	Size: 300 x 265 mm (b x h)	1. AC: 21.10.2019 bka	5. AC:
	MN 63207 MC 0000	2. AC: 24.10.2019 bka	6. AC:
	AM 33841 MZ/AZ 0000	3. AC:	7. AC:
	Typestyle: MyriadPro	4. AC:	8. AC:
Belotero Lips Contour GA Russland	Black Stanzfarbe	Type size: 6 pt	

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GB INSTRUCTIONS FOR USE FOR BELOTERO® LIPS CONTOUR

Belotero line viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in the following versions:

- 1) Belotero SOFT Lidocaine, consisting of:
 - implantable gel, in a 1.0 ml syringe - 1 pcs;
 - sterile needle 30G ½ (0,3 x 13) - 2 pcs;
 - label - 2 pcs;
 - Instructions for use – 1 pc.
- 2) Belotero BALANCE: Belotero BALANCE Lidocaine; Belotero LIPS Contour, consisting of:
 - implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc;
 - sterile needle 27G ½ (0,4 x 13) and 30G ½ (0,3 x 13) - 1 pc. each;
 - label - 2 pcs;
 - Instructions for use – 1 pc.
- 3) Belotero INTENSE: Belotero INTENSE Lidocaine; Belotero LIPS Shape, consisting of:
 - implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc;
 - sterile needle 27G ½ (0,4 x 13) - 2 pcs;
 - label - 2 pcs;
 - Instructions for use – 1 pc.

Intended Use
Belotero line in the following versions: Belotero SOFT Lidocaine, Belotero BALANCE Lidocaine, Belotero LIPS Contour, Belotero INTENSE Lidocaine, Belotero LIPS Shape is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars, filling of facial wrinkles and folds as well as for lip enhancement. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort. For use by medical practitioner. For successful treatment, it is essential that the practitioner has received a special training on the injection techniques for soft tissue augmentation and for increasing the volume of skin tissues. A good knowledge of the anatomy and physiology of the site to be treated is required.

Description
Belotero® Lips Contour is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. Belotero® Lips Contour contains 0.3% of lidocaine hydrochloride.

Presentation
Belotero® Lips Contour is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe,

two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

Composition
Cross-linked sodium hyaluronate 22.5 mg/ml
Lidocaine hydrochloride: 3.0 mg/ml
Phosphate buffer pH 7 q.s.: gel volume
The volume of the gel in each syringe is as stated on the external box.

Indications
Belotero® Lips Contour is an injectable resorbable implant indicated to correct facial atrophic scars and for lip enhancement (e.g. lip contour, perioral fine lines, moderate anti-commissures). The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

Posology and administration method
Belotero® Lips Contour is designed to be injected into the superficial to mid-dermis and the mucous membrane of the lips by a legally approved practitioner. Belotero® Lips Contour can be used for all skin types.

For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for filling wrinkles and folds and for lip enhancement. A good knowledge of the anatomy and physiology of the site to be treated is required.

The treatment must be carried out under appropriate aseptic conditions. Belotero® Lips Contour must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.

It is recommended to use one of the supplied needles. To ensure optimal use of Belotero® Lips Contour, it is recommended to assemble the needle according to the diagrams below. The practitioner should note that the use of 30G½ needles requires greater force to inject the gel and to improve patient comfort.

Belotero® Lips Contour can be used in combination with other Belotero® products during the same session. Instructions for use of each product should be followed.

If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.

Belotero® Lips Contour should be injected slowly. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.

Gently massage the treated area after the injection to distribute the product uniformly.

Contra-indications
Belotero® Lips Contour is contra-indicated:
• In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride or to amide-type local anesthetics.

• In pregnant and breast-feeding women,
• In young patients under 18 years old,
• In patients presenting a general infection.

Do not inject Belotero® Lips Contour into blood vessels. Do not inject Belotero® Lips Contour into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...). Do not inject Belotero® Lips Contour into an area previously treated with a permanent dermal filler.

Precautions for use
Before treatment the patient must be informed about the device, its contra-indications and possible side effects.

In the absence of available clinical data on tolerance and efficacy of the injection of Belotero® Lips Contour in patients with antecedents or with an active autoimmune disease, or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject Belotero® Lips Contour on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject Belotero® Lips Contour in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or cheloids. Belotero® Lips Contour can be injected in the glabellar area only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area could potentially lead to local vascular occlusion, ischemia and necrosis.

Belotero® Lips Contour can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.

Belotero® Lips Contour must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, Belotero® Lips Contour must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of Belotero® Lips Contour with the above-mentioned treatments.

No clinical data is available on the injection of Belotero® Lips Contour into an area already treated with another filling product.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject Belotero® Lips Contour on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.

For normal healthy adults, it is recommended that the maximum total dose

of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.

When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents. Do not transfer Belotero® Lips Contour into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe. Belotero® Lips Contour must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.

Discard the syringe and the remaining product after use. Do not re-sterilize and do not reuse due to the associated risks including infection.

The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

Incompatibilities
Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that Belotero® Lips Contour does not come into contact with such substances. There is no known interaction with other local or loco-regional anesthetics.

Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain sometimes accompanied by itching in the treated area. These reactions may last for a week.

- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Assembly of needle to syringe
For optimal use of Belotero® Lips Contour, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the **syringe** and the **Luer-lock adaptor** between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and unscrew it.
3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

Storage
Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

Disposal
Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with medical procedures adopted in medical institutions and local, state or federal regulations.

Requirements for environmental protection
The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

References
Updated documentation may be available from ANTEIS SA, Switzerland.

Terms of warranty
Anteis SA, 18, Chemin des Aulx/1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero line viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage, and use in accordance with the instructions for use for two years.

	Approved for printing	(Variation) Procedure No.	Identical to approved texts	Compliant with national requirements	Double checked	Date, Signature
Research and Development	<input type="checkbox"/>					
Regulatory Affairs	<input type="checkbox"/>		Yes <input type="checkbox"/>	Yes <input type="checkbox"/>	Yes <input type="checkbox"/>	
Information Officer	<input type="checkbox"/>					
Legal Department	<input type="checkbox"/>					